

FIR sauna blanket

Operation Manual



Instruction

Product feature and function

The effect of infrared slimming blanket:

- Safety-Water proof soft material with flexible tensile capabilities even in various sleeping position.
- Safety-The Heating wire is design to withstand any stretches or bends.
- Outside are top grade PU material , inside are waterproof PVC .
- Infrared dry heating.
- Two zone-temperature control.
- There are 4 pcs Thermostat which can automatically cut off for over heating.
- Throught Controller to adjust the temperature from 25-75 celsius. The time can set up from 15-60 min.
- Two zone slimming blanket are make the front and back of body to loss weight and exhause detoxin.
- The pressure on our spine is halved when you are lying down.
- With the sauna blanket ,your muscles will be loosened and you whole body relaxes.
- Your blood circulates better and the strain on your heart reduced. Our capillaries are expanded while oxygen is flowing throughout our bodies. High quality thick PU soft material.

2.product Function.


- Keep beauty-Promote blood circulation,accelerate metabolism.
Strength muscle nutrient, it also active cell and eliminate heavy ions and toxin.
- Weight loss-accelerate metabolism.so as to make the cell of fat are break down quickly.
- Body shaping –It can let skin to absorb the nutrient and inside the skin.
Make the skin smooth and beauty.
- Keep health-Improve Insomnia , constipation, endocrine dyscrasia and arthritis.Improve human immunity.

Caution

- Old people should be warded by the adults.
- Please stop using it ,when you feel bad or any other discomfortable.
- You'd better be far away from fire or other heat producer when using it.
- Please turn off and unplug the product after ysing it.

Warning! The person who has the symptom as following, please get the doctor's permission before using this product.

- The person who has product skin diseases, such as scald and reddening.
- The person who has splanchnic diseases.
- The person who has heart defibrillators or other implantable devices.
- The person who has heart disease.
- The person who has skeletal fragility or damage in his bone, especially the bone of back.
- The person who has maliqnancy.
- The person who does not feel well or is being treated.
- The person who just has a surgery.
- The person who has osteoporosis.
- The person whose blood pressure is not normal.
- The pregnant women.
- The person who has a headache.



1. Insert the cable plug of the heating blanket into the output jack in front of the main control box and fasten it. (Note: make sure power off before insertion)



2. Insert one end of the double-headed power cord into the power cord socket behind the main control box.

3. Power on:

A. Push the master power switch behind the main control box. (The device will give out auto-test sounds.)

B. The device begins to work by push the POWER key.



4. Working temperature setting:

Push the  and  keys below NUMBER A and NUMBER B. The max temperature could be set between 25° and 75°.

5. Upper/lower heating outlet selection:

Push the SELECT key, and there will be red dots at the lower right of the digitron above NUMBER A and NUMBER B, indicating that the corresponding outlet is set to output.

6. Working time setting:

Push the  and  keys above TIME. The maximum time is 60 minutes, and the minimum time is 20 minutes.

7. Heating

Push the START/STOP key, odd times for running, even times for stop. (During output, the red dots at the lower right of the digitron above NUMBER A and NUMBER B will turn on, indicating that there is output from the corresponding outlet.)